


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Methylresorcinol - technical grade, 90%

Product Number: M80407
CAS Number: 608-25-3
MDL: MFCD00002271
Formula: C₇H₈O₂
Formula Weight: 124.14 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

White to Tan to Brownish Orange

Appearance (Form)

Powder or Crystals

Infrared Spectrum

Conforms to Structure

Purity (GC)

≥ 88.5 %

Solubility (Turbidity)

Clear to Slightly Hazy

Solubility (Color)

Colorless to Orange

10% in Dioxane

Specification: PRD.1.ZQ5.10000028657

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.