


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Polypropylene - Isotactic, average M_w ~12,000, average M_n ~5,000

Product Number: 428116
CAS Number: 9003-07-0
MDL: MFCD00084447
Formula: C₃H₆



TEST	Specification
Appearance (Color)	Conforms to Requirements
Colorless to White	
Appearance (Form)	Pellets
Infrared Spectrum	Conforms to Structure
Miscellaneous Assay	160.0 - 170.0 °C
Softening Point, Ring & Ball	
Viscosity	250 - 950 cps
at 190 Deg C	

Specification: PRD.2.ZQ5.10000018719

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.