

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Strontium titanate - powder, 99%

Product Number: 396141
CAS Number: 12060-59-2
MDL: MFCD00049554
Formula: O3SrTi
Formula Weight: 183.49 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

Off-White to Beige or Grey

Appearance (Form)

Powder

X-Ray Diffraction

Conforms to Structure

Loss on Ignition

≤ 1.00 %

1 Hour at 1000 Degrees Celsius

Specification: PRD.0.ZQ5.10000050044

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.