

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
ColorBurst™ Electrophoresis Marker – mol wt 8,000–220,000 Da

Product Number: C1992

MDL: MFCD03456450

Storage Temperature: -20 °C

TEST	Specification
Appearance (Color)	Dark Purple
Appearance (Form)	Liquid
SDS Electrophoresis	Pass
Single band observed for each of the 8 proteins when 10 ul is loaded per lane.	
Western Blot	Pass
Eight major protein bands observed when 5 ul and 3 ul are transferred to nitrocellulose membrane from a tris-glycine gel	
Recommended Retest Period	Complete
1 year	

Specification: PRD.1.ZQ5.10000020534

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.