


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(Hydroxypropyl)methyl cellulose - average $M_n \sim 86,000$

Product Number: **423203**
CAS Number: 9004-65-3
MDL: MFCD00131360



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Loss on Drying	≤ 5.0 %
Miscellaneous Assay Methoxyl	28.0 - 30.0 %
Miscellaneous Assay Hydroxypropoxyl	7.0 - 12.0 %
Viscosity c = 2%; Water	3000 - 5600 cps

Specification: PRD.1.ZQ5.10000005757

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.