


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Glyceryl trioleate - $\geq 99\%$

Product Number: T7140
CAS Number: 122-32-7
MDL: MFCD00137563
Formula: C57H104O6
Formula Weight: 885.43 g/mol
Storage Temperature: -20 °C



TEST	Specification
------	---------------

Appearance (Color)	Colorless to Very Faint Yellow
Appearance (Turbidity)	Clear
Appearance (Form)	Oil
Infrared spectrum	Conforms to Structure
Purity (TLC)	$\geq 99\%$
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
100 mg/ml, CHCl ₃	
Purity (GC)	$\geq 99\%$
Recommended Retest Period	-----
3 Years	

Specification: PRD.0.ZQ5.10000043476

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.