

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

SIGMAFAST™ Protease Inhibitor Cocktail Tablets, EDTA-Free – for use in purification of Histidine-tagged proteins

Product Number: **S8830**

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	White to Off-White
Appearance (Form)	Round Tablet
Solubility (Turbidity)	Clear to Slightly Hazy
Solubility (Color)	Colorless
Solubility (Form)	Solution
Suitability	Suitable

Protease Inhibition: A 10X concentrate solution of Inhibitor will inhibit a 1.62 mg/ml suspension of P1625 &gt; or = 95%.

Specification: PRD.4.ZQ5.10000017345

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.