


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
N-(Phosphonomethyl)glycine - plant cell culture tested, BioReagent

Product Number: P9556
CAS Number: 1071-83-6
MDL: MFCD00055350
Formula: C₃H₈NO₅P
Formula Weight: 169.07 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear to Slightly Hazy
50 mg/mL, 1 N NaOH	
Infrared spectrum	Conforms to Structure
Phosphorus (P)	17.4 - 19.2 %
Plant Cell Culture Test	Pass

Specification: PRD.0.ZQ5.1000063311

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.