

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Grubbs Catalyst, 1st Generation - 97%

Product Number: 579726
CAS Number: 172222-30-9
MDL: MFCD01090946
Formula: C₄₃H₇₂Cl₂P₂Ru
Formula Weight: 822.96 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)

Conforms to Requirements

Purple or Violet

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Phosphorus NMR Spectra

Conforms to Structure

Carbon

60.6 - 64.9 %

Specification: PRD.1.ZQ5.10000005813

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.