

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification


Product Name:
2-Methyl-2-butanol - anhydrous, ≥99%

Product Number: 721123

CAS Number: 75-85-4

Formula: C₅H₁₂O

Formula Weight: 88.15 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 99.0 %
Water (by Karl Fischer)	≤ 0.003 %
Water specification for all unit sizes greater than 100 mL is 0.003%.	
Water specification for the 100 mL unit size is 0.005%	
Residue on Evaporation	≤ 0.0005 %

Specification Date : 05/02/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.