


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Lithium diisopropylamide - 97%

Product Number: 246611
CAS Number: 4111-54-0
MDL: MFCD00064449
Formula: C₆H₁₄LiN
Formula Weight: 107.12 g/mol



TEST

Specification

Appearance (Color)	White to Faint Beige
Appearance (Form)	Powder or Powder with Chunk(s)
Titration Sec-Butanol	≥ 96.5 %

Specification: PRD.0.ZQ5.10000041638

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.