

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

N-[3-(Trimethoxysilyl)propyl]-N'-(4-vinylbenzyl)ethylenediamine hydrochloride solution - 40 wt. % in methanol

Product Number:

414158

CAS Number:

33401-49-9

MDL:


MFCD00054284

Formula:

C₁₇H₃₀N₂O₃Si.ClH

Formula Weight:

374.99 AMU

**TEST****Specification**

Appearance (Color)

Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Refractive index at 25 °C

1.3900 - 1.4050

Chloride Content (in IPA)
in IPA

3.40 - 3.80 %

Specification: PRD.0.ZQ5.10000017988

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.