


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Trimethylsilyl 2,2-difluoro-2-(fluorosulfonyl)acetate

Product Number: 531421
CAS Number: 120801-75-4
MDL: MFCD02093343
Formula: C₅H₉F₃O₄SSi
Formula Weight: 250.27 g/mol



TEST	Specification
Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Proton NMR spectra	Conforms to Structure
Purity (GC)	≥ 90.0 %

Specification: PRD.1.ZQ5.10000046033

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.