


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Furanylboronic acid - ≥95.0%

Product Number: 464910
CAS Number: 13331-23-2
MDL: MFCD00799544
Formula: C₄H₅BO₃
Formula Weight: 111.89 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)

Conforms to Requirements

White to Grey to Brown

Appearance (Form)

Conforms to Requirements

Powder, Crystals, Crystalline Powder and/or Chunks

Proton NMR Spectrum

Conforms to Structure

Purity (Titration by NaOH)

≥ 95.0 %

Recommended Retest Period

6 Months

Specification: PRD.1.ZQ5.10000037970

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.