

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Allyl bromide - reagent grade, 97%, contains ≤ 1000 ppm propylene oxide as stabilizer

Product Number: 337528
CAS Number: 106-95-6
MDL: MFCD00000244
Formula: C₃H₅Br
Formula Weight: 120.98 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 96.5 %
Stabilizer	≤ 1000 ppm
Propylene Oxide	

Specification: PRD.0.ZQ5.10000024288

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.