


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: M5319
CAS Number: 519-02-8
MDL: MFCD00210527
Formula: C₁₅H₂₄N₂O
Formula Weight: 248.36 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)	White to Light Yellow
Appearance (Form)	Powder
Proton NMR spectrum	Conforms to Structure
Titration by H ₂ SO ₄	≥ 97 %
Specific Rotation	31 - 36 °
C= 1 in EtOH at 20 deg C	
Solubility (Color)	Colorless to Light Yellow
Solubility (Turbidity)	Clear
50 mg/mL, H ₂ O	

Specification: PRD.0.ZQ5.10000050178

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.