

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-B1 Bradykinin Receptor antibody produced in rabbit – affinity isolated antibody, buffered aqueous solution

Product Number: **B5560**

MDL: MFCD06798426

TEST	Specification
Formulation Solution in phosphate buffered saline containing 0.1% sodium azide.	Conforms
Concentration 1 mg/mL	Conforms
Working Dilution by Immunohistology 5 to 9 ug/mL, using human brain neurons	Conforms

Specification: PRD.0.ZQ5.10000029754

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.