


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
2-Phenylpiperazine - 96%

**Product Number:** 638641  
**CAS Number:** 5271-26-1  
**MDL:** MFCD01871362  
**Formula:** C<sub>10</sub>H<sub>14</sub>N<sub>2</sub>  
**Formula Weight:** 162.23 g/mol



TEST	Specification
------	---------------

Appearance (Color)

Conforms to Requirements

White to Orange to Tan

Appearance (Form)

Conforms to Requirements

Powder, Crystals and/or Chunks

Infrared spectrum

Conforms to Structure

Purity (Titration by HClO<sub>4</sub>)

95.5 - 104.5 %

Purity (GC)

≥ 95.5 %

Specification: PRD.0.ZQ5.10000030731

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.