


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Naphthol - 99%

Product Number: 185507
CAS Number: 135-19-3
MDL: MFCD00004067
Formula: C₁₀H₈O
Formula Weight: 144.17 g/mol



TEST

Specification

Appearance (Color)	Off White to Beige and Faint Brown to Light Brown
Appearance (Form)	Powder or Flakes
Melting Point	121 - 125 °C
Infrared spectrum	Conforms to Structure
GC (area %)	≥ 98.5 %

Specification: PRD.1.ZQ5.10000022203

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.