


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Benzoin methyl ether - 96%

**Product Number:** B8703  
**CAS Number:** 3524-62-7  
**MDL:** MFCD00008492  
**Formula:** C<sub>15</sub>H<sub>14</sub>O<sub>2</sub>  
**Formula Weight:** 226.27 g/mol



### TEST

### Specification

Appearance (Color)  
Appearance (Form)  
Powder, Crystalline Powder, Crystals, Granules and/or Chunks  
Infrared spectrum  
Purity (GC)

White to Light Yellow  
Conforms to Requirements  
  
Conforms to Structure  
≥ 95.5 %

Specification Date : 06/21/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.