


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Guaiacol - oxidation indicator

Product Number: **G5502**
CAS Number: 90-05-1
MDL: MFCD00002185
Formula: C₇H₈O₂
Formula Weight: 124.14 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Colorless to Light Yellow and Pink or White to Light Yellow and Pink
Appearance (Form)	Liquid or Solid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.0 %

Specification: PRD.0.ZQ5.10000026654

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.