

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
EX-CELL® Advanced CHO Feed 1 (with glucose)

Product Number: 24367C

Storage Temperature: 2 - 8 °C

TEST	Specification
Origin	Non-Animal Source
Bioburden	≤ 100 CFU/g
Osmolality w/o NaHCO <sub>3</sub>	675 - 775 mOs/kg
Solubility Soluble at 1X	Pass
Endotoxin Level	≤ 10 EU/ml
Site of Manufacturing	Country of Origin
Appearance (Color)	Off-White to Tan
Appearance (Form)	Free Flowing Powder

Specification: PRD.3.ZQ5.10000057989

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.