


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ethyltriphenylphosphonium bromide - 99%

Product Number: E50604
CAS Number: 1530-32-1
MDL: MFCD00011838
Formula: C₂₀H₂₀BrP
Formula Weight: 371.25 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Conforms to Requirements
Powder or Crystals	
Proton NMR spectra	Conforms to Structure
Titration by AgNO ₃	98.5 - 101.5 %

Specification: PRD.1.ZQ5.10000033078

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.