

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Poly(ethyleneimine) solution – analytical standard, 50 % (w/v) in H₂O

Product Number:

P3143

MDL:

MFCD00084427



TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid
Nitrogen	12.0 - 21.0 %

Specification: PRD.0.ZQ5.10000100965

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.