


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Ethyl-2-phenylmalonamide monohydrate - 99%

Product Number: 195022
CAS Number: 80866-90-6
MDL: MFCD00149557
Formula: C₁₁H₁₄N₂O₂ · H₂O
Formula Weight: 224.26 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Carbon	58.0 - 59.8 %
Nitrogen	12.2 - 12.8 %
Purity (GC)	≥ 98.5 %

Specification: PRD.0.ZQ5.10000043504

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.