

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

4-Imidazolecarboxaldehyde 98%

Product Number:

456128

CAS Number:

3034-50-2

MDL:


MFCD00173726

Formula:

C₄H₄N₂O

Formula Weight:

96.09 g/mol



TEST

Specification

Appearance (Color)

White to Yellow

Appearance (Form)

Conforms to Requirements

Powder or Crystalline Powder

Infrared spectrum

Conforms to Structure

Purity (Titration by HCL04)

97.5- 102.5 %

Purity (GC)

≥97.5 %

Remarks:

Specification Date : 10/25/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.