

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1-(OXAZOLO[4,5-B]PYRIDIN-2-YL)-6-PHENYL

Product Number:

ALD00032

CAS Number:

288862-83-9

Formula:

C₁₈H₁₈N₂O₂

Formula Weight:

294.35 g/mol

Storage Temperature:

-20 °C

TEST

Specification

Appearance (Color)

White to Yellow and Faint Orange to

Appearance (Form)

Orange and Faint Beige to Beige and Faint Brown to Brown

Proton NMR Spectrum

Powder or Solid or Crystals

Conforms to Structure

Specification: PRD.0.ZQ5.10000053423

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.