

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Di(propylene glycol) butyl ether, mixture of isomers ≥98.5%

Product Number:

484237

CAS Number:

29911-28-2

MDL:

MFCD00192112

Formula:

C₁₀H₂₂O₃

Formula Weight:

190.28 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥98.5 %
Water (by Karl Fischer)	≤0.30 %
Color Test	≤20 APHA
Supplier information Product of DOW Chemical	Confirmed

Remarks:

Specification Date : 06/22/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.