

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: 244953
CAS Number: 123334-22-5
MDL: MFCD00150524
Formula: C₂H₄Cl₃KPt · xH₂O
Formula Weight: 368.59 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
------	---------------

Appearance (Color)	Yellow to Yellow-Orange
Appearance (Form)	Conforms to Requirements
Powder, Fibers or Crystals	
Infrared spectrum	Conforms to Structure
ICP Major Analysis	Confirmed
Confirms Potassium and Platinum Components	
Assay	≤ 51.7 %
Platinum	

Specification: PRD.0.ZQ5.10000030217

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.