

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1-Phenyl-1-cyclopropanecarboxylic acid - 97%

Product Number:

160075

CAS Number:

6120-95-2

MDL:


MFCD00001288

Formula:

C₁₀H₁₀O₂

Formula Weight:

162.19 g/mol

**TEST****Specification**

Appearance (Color)

White to Tan

Appearance (Form)

Powder or Crystals

Infrared Spectrum

Conforms to Structure

GC (area %)

≥ 96.5 %

Specification: PRD.0.ZQ5.10000042551

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.