

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1-Phenylvinylboronic acid pinacol ester - 96%

Product Number:

659193

CAS Number:

143825-84-7

MDL:


MFCD06659923

Formula:

C₁₄H₁₉BO₂

Formula Weight:

230.11 g/mol



TEST

Specification

Appearance (Color)

White to Off White

Appearance (Form)

Crystal to Powder

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 95.5 %

Proton NMR spectra

Conforms to Structure

Specification: PRD.0.ZQ5.10000017272

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.