

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tris-Glycine Buffer 10x Concentrate

Product Number: T4904

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Turbidity)

Clear

Appearance (Color)

Colorless

Appearance (Form)

Liquid

pH

8.1 - 8.5

Protease Detection

None Detected

Suitability

Suitable

Suitable for use in gel electrophoresis and blotting.

Specification: PRD.1.ZQ5.10000007615

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.