


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Poly(tetrafluoroethylene) – powder (free-flowing), ≤12 μm particle size

Product Number: 430943
CAS Number: 9002-84-0
MDL: MFCD00084455
Formula: (CF₂CF₂)_n



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Particle Size
 micron

White to Off-White
Powder
Conforms to Structure
≤ 12

Specification: PRD.0.ZQ5.10000020898

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.