

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Wintergreen oil - China origin

Product Number: **W311308**
CAS Number: 90045-28-6
MDL: MFCD00002214



| TEST | Specification |
|----------------------------|--------------------------|
| Appearance (Color) | Colorless, Yellow or Red |
| Appearance (Form) | Liquid |
| Refractive index at 20 ° C | 1.535 - 1.538 |
| Infrared spectrum | Conforms to Structure |
| Purity (GC) | ≥ 98.0 % |
| Specific Gravity | 1.176 - 1.182 |
| Solubility (Turbidity) | Pass |
| 1ml/7ml, 70% Ethanol | |
| Acid Value | ≤ 1.0 ml |
| Arsenic (As) | ≤ 3 ppm |
| Cadmium (Cd) | ≤ 1 ppm |
| Mercury (Hg) | ≤ 1 ppm |
| Lead (Pb) | ≤ 10 ppm |
| Expiration Date Period | ----- |
| 5 Years | |

Specification: PRD.1.ZQ5.10000008098

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.