

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,5-Naphthalenedisulfonic acid disodium salt hydrate - 95%

Product Number: **250899**
 CAS Number: 207569-02-6
 MDL: MFCD00064178
 Formula: C₁₀H₆Na₂O₆S₂ · xH₂O
 Formula Weight: 332.26 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃ % NaCl	≤ 2.0 %
Water (by Karl Fischer) Typically 3-6% H ₂ O	≤ 10.00 %
Purity (TLC)	≥ 95 %

Specification: PRD.0.ZQ5.10000020020

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.