


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Sulfophthalic acid trisodium salt - technical grade

Product Number: **292362**
CAS Number: 3325-08-4
MDL: MFCD00064384
Formula: C₈H₃Na₃O₇S
Formula Weight: 312.14 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum

White to Yellow
Powder or Crystal or Chunk(s)
Conforms to Structure

Specification: PRD.1.ZQ5.10000038133

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.