

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Pyridinium dichromate on silica gel - 200-400 mesh, extent of labeling: 20 wt. % silica gel loading

Product Number:

569828

MDL:

MFCD04041078



Storage Temperature:

2 - 8 °C

TEST**Specification**

Appearance (Color)

Faint Yellow to Dark Yellow and Faint Orange to Dark Orange

Appearance (Form)

Powder

Loss on Drying

≤ 4.0 %

110 Degrees Celsius

Loading, Molecular

18.0 - 22.0 %

% Product

Size

Confirmed

200-400 Mesh

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.