


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
trans-2-Hydroxy-β-nitrostyrene - 97%

Product Number: 639680
CAS Number: 3156-43-2
MDL: MFCD01204157
Formula: C₈H₇NO₃
Formula Weight: 165.15 g/mol



TEST	Specification
Appearance (Color)	White to Dark Yellow and Faint Beige to Beige to Orange
Appearance (Form)	Conforms to Requirements
Powder, Crystals and/or Chunks	
Infrared spectrum	Conforms to Structure
Carbon	56.1 - 60.2 %
Nitrogen	8.2 - 8.8 %
GC (area %)	≥ 96.5 %

Specification: PRD.3.ZQ5.10000033157

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.