


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Poly(Bisphenol A carbonate) – analytical standard

Product Number: 181641
CAS Number: 25037-45-0
MDL: MFCD00084476
Formula: C₁₅H₁₆O₂



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Pellets
Infrared spectrum	Conforms to Structure
Molecular Weight by Light Scattering	
Typically 26,000 to 30,000	
Molecular Weight by GPC	
Typically 26,000 - 30,000	
Molecular Number	
Typically 16,000 to 18,000	
Molecular Wt/Molecular Nm	
Typically 1.44 to 1.88	

Specification: PRD.1.ZQ5.10000026470

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.