

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium nitrite - ReagentPlus[®], ≥99.0%

Product Number: S2252
CAS Number: 7632-00-0
MDL: MFCD00011118
Formula: NNaO₂
Formula Weight: 69.00 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder or Crystals
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless to Faint Yellow
Purity	≥ 99.0 %
Recommended Retest Period	-----
2 Years	

Specification: PRD.0.ZQ5.10000023761

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.