


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Isocytosine - ≥99%

Product Number: **I2127**
CAS Number: 108-53-2
MDL: MFCD00057557
Formula: C₄H₅N₃O
Formula Weight: 111.10 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity) 50 mg/mL in HOAc	Clear
EMM	6.40 - 6.70
Lambda max 273 nm in 0.1 M NaOH	
Purity (TLC)	≥ 99 %

Specification: PRD.0.ZQ5.10000061723

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.