

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Ammonium iron(II) sulfate hexahydrate - ReagentPlus®, ≥98%

**Product Number:** F3754  
**CAS Number:** 7783-85-9  
**MDL:** MFCD00150530  
**Formula:** H8FeN2O8S2 · 6H2O  
**Formula Weight:** 392.14 g/mol

 $(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$ 

TEST	Specification
------	---------------

Appearance (Color)	Light Green to Very Dark Green and Faint Blue to Dark Blue and Blue-Green
Appearance (Form)	Powder or Crystals
Solubility (Color)	Faint Blue to Dark Blue and Faint Green to Very Dark Green and Very Faint Blue-Green to Blue-Green
Solubility (Turbidity)	Clear to Very Slightly Hazy
10G plus 99 mL H2O and 1 mL H2SO4	
Titration by KMNO4	≥ 98 %

Specification: PRD.3.ZQ5.10000021724

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.