


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
TBS-DHG Catalyst - ≥95%

Product Number: **901235**
CAS Number: 212120-61-1
Formula: C₁₂H₂₆O₄Si
Formula Weight: 262.42 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals or Solid or Chunks
¹ H NMR Spectrum	Conforms to Structure
Carbon Content	52.2 - 57.7 %
Specific Rotation 0.1%, CHCl ₃	-25.0 - -15.0 °
Purity (GC)	≥ 95.0 %

Specification: PRD.0.ZQ5.10000088703

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.