


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
6-Methoxysalicylic acid - 98%

Product Number: 464139
CAS Number: 3147-64-6
MDL: MFCD00674090
Formula: C₈H₈O₄
Formula Weight: 168.15 g/mol



TEST

Specification

Appearance (Color)	White
Appearance (Form)	Powder or Solid or Crystals
Infrared spectrum	Conforms to Structure
Purity (Titration by NaOH)	97.5 - 102.5 %
Solubility (Turbidity)	Clear
c = 25mg/ml, Methanol	
Solubility (Color)	Colorless
Purity (TLC)	≥ 98 %

Specification: PRD.0.ZQ5.10000038112

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.