

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Product Number: 791334



TEST	Specification
Appearance (Color)	Colorless or White
Appearance (Form)	Liquid
ICP Major Analysis Confirms Silicon Component	Confirmed
pH	1 - 3
Solid Content	24.0 - 30.0 %
Particle Size (Dv 50)	≤ 30 nm
Recommended Retest Period 1 Year	-----

Specification: PRD.0.ZQ5.10000097278

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.