

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Molybdenum disulfide ink for spin/spray coating

Product Number: **901867**

TEST	Specification
Appearance (Color)	Green and Green-Grey
Appearance (Form)	Suspension
Solid Content	1.8 - 3.0 %
ICP Major Analysis	Confirmed
Confirms Mo and S Components	
Recommended Retest Period	-----
9 Months	

Specification: PRD.0.ZQ5.10000095922

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.