

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)


## Product Specification

Product Name:

Cortisone-2,3,4-<sup>13</sup>C<sub>3</sub> solution - 100 µg/mL in methanol, 98 atom % <sup>13</sup>C, 97% (CP)

Product Number:

803154



Storage Temperature:

-20 °C

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Appearance (Clarity)	Clear
Purity (HPLC)	≥ 97 %
Proton NMR Spectrum	Conforms
Carbon NMR Spectrum	Conforms
Mass Spectrometry-13C	≥ 98 %
Concentration	90 - 110 µg/ml

Specification: PRD.0.ZQ5.10000058807

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.