

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Silver trifluoromethanesulfonate - ≥99%

Product Number:**176435**

CAS Number:

2923-28-6

MDL:


MFCD00013226

Formula:

CAgF3O3S

Formula Weight:

256.94 g/mol

**TEST****Specification**

Appearance (Color)

White

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Titration by KSCN

41.5 - 42.4 %

% Silver

Solubility (Turbidity)

Clear to Slightly Hazy

50mg/ml CH₃CN

Solubility (Color)

Colorless

Specification: PRD.1.ZQ5.10000023680

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.