

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:


Potassium cyclopropyltrifluoroborate

Product Number:

662984

CAS Number:

1065010-87-8



Formula:

C3H5BF3K

Formula Weight:

147.98 g/mol

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Proton NMR spectra	Conforms to Structure
Boron HNMR Spectra	Conforms to Structure
Fluorine NMR Spectra	Conforms to Structure
Purity Titration NAOH by HCL	88.0- 112.0 %
Carbon	21.43- 27.27 %

Remarks:

Specification Date : 11/11/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.