

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Cyanophenylboronic acid - $\geq 95.0\%$

Product Number: 513016
CAS Number: 150255-96-2
MDL: MFCD01318967
Formula: C₇H₆BNO₂
Formula Weight: 146.94 g/mol

TEST	Specification
Appearance (Color)	Conforms to Requirements
White to Off-White to Light Yellow to Light Orange	
Appearance (Form)	Powder
Proton NMR Spectrum	Conforms to Structure
Purity (Titration by NaOH)	$\geq 95.0\%$

Specification: PRD.1.ZQ5.10000015636

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.